

## Placenta Ribonuclease Inhibitor

<b>Catalog No.:</b>	<b>Z5040001</b>
<b>Lot No.:</b>	
<b>Source:</b>	Human Placenta
<b>Size:</b>	5,000 units
<b>Concentration:</b>	40 u/μl
<b>Storage Buffer:</b>	20 mM HEPES-KOH (pH 7.5), 50 mM KCl, 5 mM DTT and 50% glycerol.
<b>Unit Definition:</b>	One unit is the amount of enzyme required to inhibit by 50% the activity of 5 ng of RNase A at 25°C (This inhibitor activity is determined by its ability to inhibit hydrolysis of cyclic 2', 3'-CMP by RNase A)
<b>Quality Control:</b>	Nuclease activity was not detected in any of the following cases, as judged from the intact gel electrophoresis: <ol style="list-style-type: none"><li>1. After incubation of 1 μg of DNA-Hind III fragments with 300 units of enzyme for 24 hours at 37°C</li><li>2. After incubation of 1 μg of pBR322 DNA with 300 units of enzyme for 24 hours at 37°C</li><li>3. Migrates as a single band of 50 kd in SDS-polyacrylamide gel electrophoresis.</li></ol>

Note: Enzyme inhibits in a wide pH range, but most strongly at pH 7-8. Requires DTT of at least 1 mM to be active.

**Storage and Handling:** -20°C

**Shelf Life:** One year from the date of purchase if stored and handled properly

## Protocol

### A. Transcription in vitro (unlabeled RNA)

The standard in vitro transcription assay below uses Ribonuclease Inhibitor at a final concentration of 1 u/μl. With appropriate modifications, this reaction can be used for in vitro transcription analysis in a variety of experimental applications.

5X transcription buffer	20 μl
DTT, 100mM	10 μl
Ribonuclease Inhibitor	100 u
ATP, GTP, CTP and UTP, 2.5 mM each	20 μl
linearized plasmid DNA, 2–5 μg in H <sub>2</sub> O or TE buffer	2 μl
RNA Polymerase; SP6, T3 or T7	50 u
<u>Nuclease-free water</u>	<u>X μl</u>
Final volume	100 μl

Incubate for 60–120 minutes at 37–40°C.

### B. Transcription in vitro (<sup>32</sup>P-labeled RNA probes)

5X transcription buffer	4 μl
DTT, 100mM	2 μl
Ribonuclease Inhibitor	20 u
ATP, GTP and UTP, 2.5 mM each	4 μl
CTP, 100 μM	2.4 μl
Linearized template DNA, 0.2–1.0mg/ml in H <sub>2</sub> O or TE buffer	1 μl
[α- <sup>32</sup> P]CTP, 50μCi at 10mCi/ml	5 μl
RNA Polymerase, SP6, T3 or T7	1 μl
<u>Nuclease-free water</u>	<u>X μl</u>
Final volume	20 μl

Incubate for 60 minutes at 37–40°C.

## III. Composition of Buffers and Solutions

### 5X transcription buffer

200 mM Tris-HCl (pH 7.5)

30 mM MgCl<sub>2</sub>

10 mM spermidine

50 mM NaCl

### 1X TE buffer

10 mM Tris-HCl (pH 8.0)

All products are for research only

1 mM EDTA