

Data Sheet

Product Name: Lymphoma Tissue Array - 39 cases of various types of nodal and extranodal lymphoma


Catalog No.: Z7020072

Lot No.: C502161

Species: Human Mouse Rat Monkey (Rh) Guinea Pig Porcine
 Bovine Hamster Dog Monkey (Cy) Rabbit Plant

Tissue Type: Normal Adult Fetal Tumor Disease Cell Line

Tissue Array Diagram:

Cores	Size	Cut	Format	QA/QC	
40	1.1 mm	4 µm	5X8	H&E, IHC anti-LCA	

Recommended applications: For research use only. RNA or protein analyses using IHC or ISH suitable for general lymphoma tissue profiling.

Description: Lymphoma array, >30 cases of various types of nodal and extranodal lymphoma.

All the tissues were from surgical resection. They were fixed in 10% neutral buffered formalin for 24 hours and processed using identical SOPs. Sections were picked onto Superfrost Plus or APES coated Superfrost slides. They all have a guaranteed six months' shelf-life at 4°C from the date of shipment. There may be >5% tissue core loss, but the core retaining rate should be >90%.

Position	Sex	Age	Site	Histology	Grade	Stage
A01	F	77	Thyroid	Lymphoma, MALT-B cell lymphoma, high grade	null	null
A02	F	24	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
A03	F	24	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
A04	M	53	Intestine, small intestine	Lymphoma, peripheral T cell lymphoma	null	null
A05	M	55	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
A06	M	54	Shoulder	Lymphoma, non-Hodgkin B cell lymphoma	null	null
A07	M	39	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
A08	M	52	Leg	Lymphoma, non-Hodgkin B cell lymphoma	null	null
B01	F	27	Lymph node	Lymphoma, diffuse large B cell lymphoma	null	null
B02	M	48	Intestine, colon	Lymphoma, non-Hodgkin lymphoma	null	null
B03	M	60	Intestine, colon	Lymphoma, non-Hodgkin B cell lymphoma	null	null
B04	M	72	Peritoneal cavity	Lymphoma, non-Hodgkin B cell lymphoma	null	null
B05	M	42	Lymph node, axillary	Lymphoma, non-Hodgkin B cell lymphoma	null	null
B06	M	7	Lymph node, neck	Lymphoma, Hodgkin lymphoma	null	null
B07	M	66	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
B08	M	33	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
C01	M	51	Lymph node, neck	Lymphoma, anaplastic large cell lymphoma	null	null
C02	F	32	Lymph node, lung	Lymphoma, non-Hodgkin B cell lymphoma	null	null
C03	M	72	Testis	Lymphoma, non-Hodgkin B cell lymphoma	null	null
C04	F	71	Lymph node	Lymphoma, non-Hodgkin B cell lymphoma	null	null
C05	M	43	Stomach	Lymphoma, MALT-B cell lymphoma, high grade	null	null
C06	F	52	Lymph node, neck	Lymphoma, diffuse large B cell lymphoma	null	null
C07	M	60	Intestine, colon	Lymphoma, non-Hodgkin B cell lymphoma	null	null

C08	M	30	Lymph node, neck	Lymphoma, Hodgkin lymphoma	null	null
D01	F	33	Stomach	Lymphoma, MALT-B cell lymphoma, high grade?	null	null
D02	F	70	Lymph node, mesentery	Lymphoma, non-Hodgkin B cell lymphoma	null	null
D03	F	43	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
D04	M	76	Tonsil	Lymphoma, non-Hodgkin B cell lymphoma	null	null
D05	M	65	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
D06	M	39	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
D07	M	53	Lymph node, neck	Lymphoma, diffuse large B cell lymphoma	null	null
D08	M	64	Lymph node, neck	Lymphoma, non-Hodgkin B cell lymphoma	null	null
E01	F	58	Lymph node	Lymphoma, non-Hodgkin B cell lymphoma	null	null
E02	M	43	Lymph node, neck	Lymphoma, mantle cell lymphoma	null	null
E03	M	45	Pelvic cavity	Lymphoma, mantle cell lymphoma	null	null
E04	M	45	Pelvic cavity	Lymphoma, mantle cell lymphoma	null	null
E05	M	45	Pelvic cavity	Lymphoma, mantle cell lymphoma	null	null
E06	F	52	Lymph node	Lymphoma, mantle cell lymphoma	null	null
E07	M	69	Lymph node, neck	Lymphoma, mantle cell lymphoma	null	null
E08	M	69	Lymph node, neck	Lymphoma, mantle cell lymphoma	null	null

Notes: Bake at 60°C for 30 minutes before use. If antigen retrieving is needed, it is always a good idea to start with a protocol with weak to mild strength.

Clinical Staging: Malignant lymphoma

Stage I - Involvement of a single lymph node region (I), or localized involvement of a single extrahepatic organ or site (I_E)

Stage II - Involvement of two or more lymph node regions on the same site of the diaphragm (II), or localized involvement of a single extralymphatic organ or site and its regional lymph node(s) with or without involvement of other lymph node regions on the same side of the diaphragm (II_E). The number of lymph node regions involved may be indicated by a subscript (e.g. II₃)

Stage III - Involvement of lymph node regions on both sides of the diaphragm (III), which may also be accompanied by localized involvement of an associated extralymphatic organ or sites (III_E), or by involvement of spleen (III_s) or both (III_E+s)

Stage IV - Disseminated (multifocal) involvement of one or more extralymphatic organs, with or without associated lymph node involvement; or isolated extralymphatic organ involvement with distant (non-regional) nodal involvement

A and B Classification

Each stage should be divided into A and B according to the absence or presence of defined general symptoms. These are:

1. Unexplained weight loss of more than 10% of the usual body weight for the 6 months prior to first attendance.
2. Unexplained fever with temperature above 38 degree
3. Night sweat

Pathological Staging

Splenectomy, liver biopsy, lymph node biopsy and marrow biopsy are mandatory for the establishment of pathological staging

TNM classification of malignant tumours, Fifth Edition (1997)

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: _____

