

Certificate of Analysis

Product Name: xTractPure™ FFPE RNA Extraction Kit

Catalog No.: Z2212003

Lot No.:

Shipping Condition: Ship at Room Temperature (dry ice for DNase-Treatment box if purchased)

Storage Condition: Upon receipt, store Proteinase K at -20°C. If purchased, also store DNase-Treatment box components at -20°C

Shelf Life: 1 year from the date of receipt under proper storage conditions

Description

Components in this kit are prepared with pure chemicals according to proprietary technology. They are designed for the extraction and purification of FFPE RNA from formalin-fixed paraffin-embedded tissue samples.

Quality Control

- This kit is tested by extracting FFPE RNA from a selected tissue block. The resulting eluted RNA is tested for yield, fragment sizes and purity and compared to known yields, fragment sizes, and purity of the same control tissue block to ensure efficacy and reproducibility.

Kit Components

1. xTractPure™ FFPE DR Extraction Kit
| | With DNase-Treatment box (Z2212099)
2. Certificate of Analysis

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY:  _____