

## Certificate of Analysis

**Product Name:** cfPure® / cfPure® MAX cfDNA Extraction Kit

**Catalog No.:**  K5011610  K5011625  K5011625MA

**Lot No.:**

**Shipping and Storage Condition:** Ship and store at Room Temperature

**Shelf Life:** 1 year from the date of receipt under proper storage conditions

### Description

Components in this kit are prepared with pure chemicals according to the proprietary technology. They are designed for the extraction and purification of cell-free DNA from human plasma. This kit provides the necessary reagents to extract and elute DNA from human plasma.

### Quality Control

- A representative kit of this lot is randomly selected and tested by extracting cfDNA from human control plasma. The resulting eluted cfDNA is tested for yield and purity and compared to known yields and purity of the same control human plasma to ensure efficacy and consistency.

### Kit Components

1. cfPure® / cfPure® MAX cfDNA Extraction Kit
2. Certificate of Analysis

**FOR IN VITRO RESEARCH USE ONLY**

**APPROVED BY:**  \_\_\_\_\_