

Certificate of Analysis

Material name:	cfPure® DNA Free Plasma
Part number:	Z2010001/Z2010002
Lot Number:	
Fill Volume	5/50 ml
Source:	Defibrinated Normal Source Plasma
Storage:	2-8°C for usage within 30 days or -20 ° for long term storage. Avoid repeated

DNA Testing Results

DeNovix Ultra High Sensitive Assay*: Concentration of DNA: < 2 pg/mL

Typical Product Characteristics

Total protein: 4-7 g/dL

Immunoglobulin G: 300-1700 mg/dL

Immunoglobulin A: 70-400 mg/dL

Immunoglobulin M: 10-172 mg/dL

Cholestrol, Total: 30-199 mg/dL

Triglycerides: 30-149 mg/dL

HDL Cholestrol: > 10 mg/dL

VLDL Cholestrol: 10-40 mg/dL

LDL Cholestrol: 20-99 mg/dL

pH: 7.2-7.6

Filtered through: 0.22 um

Caution

All original donor units were collected at FDA approved donor centers located in the United States. The content of this product has been tested and found negative for HBs antigen, antibodies to HIV-1/2, HCV and Syphilis. It has been tested and found non-reactive for nucleic acids of HIV, HBV and HCV virus. All testing performed with kits approved by the FDA. No detectable amount of DNA/RNA was found from the plasma including those for HTLV I&II and Zika. However, full safety precautions should be taken as in the handling of any potentially infectious body fluid. The product should be handled only by persons experienced in laboratory techniques

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY _____

