

Certificate of Analysis

Product Name: uPure Urine DNA Extraction Kit

Catalog No.: K5011196

Lot No.:

Shipping Condition: Ship at Room Temperature

Storage Condition: Upon receipt, store Proteinase K at -20°C.

Shelf Life: 1 year from the date of receipt under proper storage conditions

Description

Components in this kit are prepared with pure chemicals according to proprietary technology. They are designed for the extraction and purification of cell-free DNA and genomic DNA from human urine.

Quality Control

- This kit is tested by extracting cell-free DNA and genomic DNA from human urine. The resulting eluted DNA is tested for yield and purity and compared to known yields and purity of the same control human urine to ensure efficacy and reproducibility.

Kit Components

1. uPure Urine DNA Extraction Kit
2. Certificate of Analysis

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY:  _____