

Certificate of Analysis

Product Name: IsHyb In Situ Hybridization (ISH) Kit

Catalog No.: K2191020/K2191050

Lot No.: _____

Shipping & Storage Condition: 2 – 8 °C

Shelf Life: 2 years from the date of receipt under proper storage conditions

Description

Biochain's IsHyb In Situ Hybridization Kit provides researchers a method doing in situ hybridization with non-radioactive probes. Compared with radioactive methods, the sensitivity of this kit is the same while the procedure is simplified. This kit is suitable for digoxigenin labeled probes combined with an alkaline phosphatase color development system. Using this kit, all disadvantages of working with radioactive materials can be avoided.

Quality Control

1. Components in this kit are free of precipitate.
2. Each lot is tested by β -actin cDNA probe to have the same result as a control lot

Components

1. 1 kit
2. Certificate of analysis
3. One complete User Manual
4. SDS

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY:  _____