

# Certificate of Analysis

**Product:** Genomic DNA Extraction Kit

**Catalog Number:** K5016005

## Description

Biochain's DNA extraction kit is a convenient tool for scientist isolate high quality genomic DNA. The isolated DNA can be used for PCR amplification template, Southern Blot analysis, SNP analysis, DNA methylation research, etc. The kit contains enough reagents for isolating genomic DNA from 5 grams of tissue and cells.

## Kit Components

1. One Genomic DNA Extraction Kit
2. Certificate of Analysis
3. One User Manual
4. MSDS

**FOR IN VITRO RESEARCH USE ONLY**

**APPROVED BY:** 