

Certificate of Analysis

Product Name: DNA Storage Solution **Lot:**

Cat: LB3401010

Shipping Condition: Room Temperature

Storage Condition: Room Temperature

Shelf Life: 1 year from the date of receipt under proper storage conditions.

Description

DNA Storage Solution is used for dissolving DNA. Unlike the general TE buffer, this solution can keep DNA concentration more stable, especially for low concentration DNA. This solution is also good for DNA for sequencing purpose, helps to avoid decreasing the cluster numbers due to repeated freeze and thaw cycles.

Quality Control

DNA dissolved in DNA Storage Solution at 50 ng/ul concentration, after 3 freeze-thaw cycles, the DNA concentration is still between 45 ng/ul and 55 ng/ul.

Components

1. 10 ml DNA Storage Solution
2. Certificate of Analysis

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: _____

