

Certificate of Analysis

Product Name: SeqPure PCR Purification Kit

Catalog No: Z1510001, Z1510002, Z1510003

Lot No: _____

Shipping Condition: Room Temperature

Storage Condition: 2-8°C

Shelf Life: 1 year from the date of receipt under proper storage conditions

Description

Components of this kit are prepared with pure chemicals according to proprietary technology. It is used for purification and analysis of PCR reaction products.

Quality Control

Components of this kit are free of precipitate. The quality and purity of DNA fragments after using SeqPure PCR purification kit were tested by a spectrophotometer and electrophoresis. $A_{260/280}$ is between 1.8 and 2.0 and $A_{260/230}$ is >2.0 . The DNA is also visualized on an agarose gel to confirm recovery of fragments 100bp or above.

Components

1. One Kit
2. User Manual
3. Certificate of Analysis
4. SDS

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY:  _____