

Certificate of Analysis

Product: Rapid Gold Immuno-assay Kit for Antibody IgG to Tuberculosis (TB)

Shipping Condition: Blue ice

Storage Condition: 2-8°C

Shelf Life: One year from the date of receipt under proper storage conditions

Description

This kit uses one step method to detect anti-TB antibodies IgG in serum or plasma. The purified recombinant TB antigen and gold labeled conjugate are coated on the the surface of nitrocellulose membrane. Add blood or serum sample on the membrane, if anti-TB antibody is present in the sample, a complex of TB-anti-TB will form, the interation of the complex and gold labeled conjugate will produce pink color on the membrane. The pink color intensity will reflect the antibody intensity in the sample. The test is fast and easy to operate.

Components

1. One Kit
2. One Certificate of Analysis
3. One User Manual

FOR IN VITRO RESEARCH USE ONLY

APPROVED BY: _____

